

# FO

# F



**FOFO** | Foshan Hongfeng Co., Ltd

## Professional Manufacturer

- Recovery boots
- Air Compression Therapy System
- Calf Wireless Mini Massage device
- Massage Gun



### Foshan Hongfeng Co.,Ltd

Address: No.4 Leqiang Road ,Leping Sanshui,Foshan,528100,GuangDong,China

Tel: 0757-88791710 ext 888

0592-2097262

Email: fofo@fofomedical.com

Mob: +86 13055506651

www.fofomedical.com



www.fofomedical.com



Star Suppliers



Dedicated to manufacturing high quality medical and health care products  
Your proprietary supplier



TO KNOW MORE  
ABOUT US



# COMPANY INTRODUCTION


With over 20 years of experience in medical and Rehabilitative products manufacturing, We have numerous customers located across all regions and continents. Being manufactured from superior materials and using modern technology, our main products cover rehabilitation, respiratory, monitoring and physical therapy, including air compression therapy for DVT, accessories like air pumps, motors and compressors. All our products are CE, FDA and ISO13485 approved and come with guarantee support. Competitive prices, prompt delivery schedules and reliable services. FOFO medical will be your ideal partner!



<b>TOP 1</b> Medical air mattress & Air compression leg massager up to Alibaba	<b>40+</b> Customers in 40+ countries worldwide
<b>\$30,000,000</b> USD in sales by 2022	<b>20 years</b> of high-quality production experience
<b>22,000</b> square meters of production area	<b>1,200,000</b> Exporting 1.2 million pieces of air mattresses in 2022

# 01

## Air Compression Therapy System



### Features and Benefits

- Compression therapy increases blood flow to specific parts of the body, accelerating its natural.
- Mechanisms for delivering oxygen, nutrients, and hormones throughout the body while simultaneously removing metabolic waste. When this process is optimized, the body can speed up muscle recovery, relieve pain and inflammation, and improve athletic performance.

### How can Compression Boots help ?

- ✔ Lessen pain sensitivity. Therefore accelerating recovery by reducing muscle soreness from pressure stimuli.
- ✔ Clear metabolites Passively. compression significantly lower blood lactate levels.
- ✔ Assist with a positive adaptive response to exercise.
- ✔ Decreases Muscle Fatigue After Exercise. Increases Flexibility and skeletal muscle-oxidative stress during recovery from heavy resistance training.
- ✔ Increased Range of Motion. enhances acute range of motion with less discomfort and time.

### Applicable people

- ① Pregnant and wish to decrease fluid build up in the legs. This is particularly problematic for ladies in their third trimester
- ② Weekend warriors and athletes who have been training hard and have sore legs
- ③ People who are on their feet all day and have tired/sore legs
- ④ Those who sit at a desk all day and have fluid build up or poor circulation in their legs
- ⑤ Runners, cyclists, swimmers or crossfitters who are trying to freshen their legs prior to an event or race

Sleeve	Description	Size
Leg	Leg cuff M	94*60cm
	Leg cuff L	100*70cm
	Leg cuff XL	110*70cm Overlapping
	Leg cuff XXXL	125*76cm OverLapping
Arm	Arm cuff M	87*51cm
	Arm cuff L	90*60cm OverLapping
Waist	Waist cuff M	120*37cm
	Waist cuff L	130*37cm OverLapping
Trousers	Short pants (4,6,8 chambers)	80*76cm
	Long pants (6,8 chambers)	142*85cm
Vest/cloth	Vest (6,8 chambers)	78*72cm
Jacket	Double arm/shoulder	133*37 (pair)



6 Massage Modes

**A** One by one, single part massage, after all chambers are treated, one cycle complete



**B** One after another, when all chambers are inflated together, one cycle complete



**C** Four times inflatable massage from bottom to top



**D** All together to inflate, and then deflate together, four times as one cycle



**E** Combined B and C mode



**F** Combined A and B mode



MODEL: **F03001**



Waist



Easy operation



Arm



Effective lymphatic edema

Manual device

Specifications:

- Device size: 24cm\*20cm\*11cm
- Massage mode: A
- Pressure range: 30~250 mmHg
- Treatment time: 10~60 mins/always on
- Power consumption: 65VA
- Each chamber can be partially skipped


A Massage Modes



MODEL: **F03001B**



Waist



Convenient one-touch sleeve



Arm



Use various cuffs according to user's area

Sequential compression devices (SCDs)

Specifications:

- Device size: 24cm\*20.5cm\*20cm
- Massage mode: therapy mode
- Pressure range: 30~250 mmHg
- Treatment time: 10~60 mins/always on
- Power consumption: 65VA
- Each chamber can be partially skipped

B Massage Modes



MODEL: **F03006-F03008**



Convenient handle design



Color Choose

Home use

Specifications:

- Device size: 24cm\*20cm\*11cm
- Massage mode: 6 modes
- Pressure range: 30~250 mmHg
- Treatment time: 10~60 mins/always on
- Power consumption: 65VA
- Each chamber can be partially skipped

# Air Compression Therapy System: Digital

## Advantage



The ability to turn off individual chambers in the cuffs, e.g. to avoid injured places

**1**



Possibility to change the parameters of the treatment while the device is operating

**2**



Automatic calculation of gradient values after setting minimum and maximum values

**3**



Available for the visually impaired, set voice prompts, customizable special keyboard

**4**



System can be updated and upgraded for continuous access to the latest features

**5**



Allows immediate termination of the treatment and deflation of air from the cuffs

**6**




Real pressure measurement during the treatment.

**7**




Adopt touchable screen, easy to operate and provide a clear display

**8**



Remote control distance >5 meters



### MODEL: F03001



Use durable domestic air motors



6 Pressure modes (Circulation&Massage)

#### Specifications:

Device size: 24cm\*20.5cm\*20cm  
 Massage mode: 6 modes  
 Pressure range: 30 ~ 250 mmHg  
 Treatment time: 10~ 60 mins  
 Power consumption: 65VA  
 Each chamber can be partially skipped  
 4,6,8 chambers available

#### Specifications:

Device size: 24cm\*20cm\*11cm  
 Massage mode: 3 modes  
 Pressure range: 30~250 mmHg  
 Treatment time: 10~60 mins  
 Power consumption: 65VA  
 Each chamber can be partially skipped



Use various cuffs according to user's area



Care for your choice of foot/calf/knee/thigh

### MODEL: F03006-F03008



### MODEL: F03001B



Set pressure and time



Skip function

#### HOT Specifications:

Device size: 24cm\*20.5cm\*20cm  
 Massage mode: 6 modes  
 Pressure range: 30 ~ 250 mmHg  
 Treatment time: 10~ 60 mins  
 Power consumption: 65VA  
 Each chamber can be partially skipped  
 4,6,8 chambers available

#### Specifications:

Device size: 31cm\*16.5cm\*11cm  
 Massage mode: 6 modes  
 Pressure range: 30 ~ 250 mmHg  
 Treatment time: 10~60 mins  
 Power consumption: 65W  
 Each chamber can be partially skipped



Effective lymphatic edema



Color option

### MODEL: F03001C



## Air Compression Therapy System: **Battery**

- \* Compact, lightweight, battery-powered devices
- \* Outdoors available, postexercise immediately treatment

### MODEL: **F03002**

**HOT** Specifications:

Device size: 24cm*20.5cm*20cm
Mode: 5 modes
Pressure range: 30 ~ 250 mmHg
Treatment time: 10~ 60 mins
Power consumption: 65VA
Each chamber can be partially skipped
Battery: 2600mA/spec



### MODEL: **F03002B**

Specifications:

Device size: 25.5cm*13.5cm*10cm
Massage mode: 5 modes
Pressure range: 30 ~250 mmHg
Treatment time: 10~60 mins
Input: DC 15V/5200mA
Each chamber can be partially skipped



### MODEL: **F03002C**

Specifications:

Device size: 27cm*13cm*10cm
Massage mode: 5 modes
Pressure range: 30~250 mmHg
Treatment time: 10~60 mins
Input: DC 15V/5200mA
Each chamber can be partially skipped



## Air Compression Therapy System: **high-end custom Bluetooth**

- \* LCD touch screen
- \* Available for the visually impaired, set voice prompts, customizable special keyboard
- \* System can be updated and upgraded for continuous access to the latest features

### MODEL: **F03003**

#### Specifications:

Model No: F03003
Specification Weight: 400g
Battery capacity: 2500mAh polymer lithium battery
Leg strap size : 600*350*2.5MM
Main Device Size: 27cm*13cm*10cm
Optional Pressure : 60/80/99mmHg
Input voltage DC5V/1500-2000 mA



Time Selection

Automatic calculation of gradient values after setting minimum and maximum values

Message mode

Each chamber can be partially skipped

On/Off

Start/Pause



## Parameter comparison

Maunal			
Model	Function	Cuff	Chamber
FO3001	1 modes	Leg,extension	4/6/8
FO3001B			
FO3006-3008			

Digital			
Model	Function	Cuff	Chamber
FO3001	6 modes	Leg,arm, waist,shorts, vest, long pants, extension	4/6/8
FO3001B			
FO3006-3008			
FO3001C			6/8

Battery			
Model	Function	Cuff	Chamber
FO3002	5 modes	Leg,arm, waist,shorts, vest, long pants, extension	4/6/8
FO3002B			
FO3002C			

WIFI			
Model	Function	Cuff	Chamber
FO3003	5 modes	Leg,arm, waist,shorts, vest, long pants, extension	4/6/8

## Calf Wireless Mini Massager

# 02

\* DONT LET FATIGUE AFFECT YOUR QUALITY OF LIFE.  
GOOD QUALITY IS IN THE DETAILS.

- MODEL: **FO3002 mini**



### Specifications :

Model No:	FO3002mini
Specification Weight:	400g
Battery capacity:	2500mAH polymer lithium battery
Leg strap size :	600*350*2.5MM
Main Device Size:	14*7*2.5mm
Optional Pressure :	60/80/99mmHg
Input voltage:	DC5V/1500-2000 mA



# 03

## Massage Gun

\* available for heating & cooling function (35-65°C)

### Full Metal Body Design

Long-handle Design  
Brushless Powerful Motor 5530  
Maxium Torque:250mN.m  
Maxium:3800RPM  
Battery Spec:22.2V-2600mAH/6PCS  
Amplitude:6MM

### Design For Strongman Over 100kg

300 Gears Touch Version  
Speed:From 400RPM to 3400RPM  
by Step 100RPM  
Support Fast Change Speed



### Packing Information

Giftbox Dimension: 380X265X90MM  
Giftbox Weight: 3.8/KG 4pcs/Carton  
Carton Dimension: 385X280X395MM  
Carton Weight: 16KG



### MODEL: HF1700



### MODEL: HF200H



### Key Parameters

Low Noise(< 50dB)  
Brushless Powerful Motor  
Maxium:3200RPM  
Battery Options  
11.1V-2600mAH/3PCS  
Amplitude:10MM

### Lady's Models

Slim Handle  
Super Silence Noise<30dB  
Six Private Heads  
Dismount battery  
Share charger with smartphones

### Packing Information

Giftbox Dimension: 285X228X92MM  
Giftbox Weight: 1.3/KG 20PCS/Carton  
Carton Dimension: 935X300X250MM  
Carton Weight: 14KG

### MODEL: HF2100H



### Key Parameters

Low Noise(40-50dB)  
Brushless Powerful Motor Maxium:3700RPM  
Battery 22.2V-2600MAH-6PCS 18650  
Amplitude:10MM,12MM,14MM for your choice

### 20 Gears-Digital LED Version Friendly software defination

Any button to start from gear 1  
Fast change speed  
Speed loop  
Gear from 1000 to 3000RPM by step 100  
Heating/Cold Head as optional

### Packing Information

Giftbox Dimension: 242X100X345MM  
Giftbox Weight: 1.8/kg 7PCS/Carton  
Carton Dimension: 750X305X268MM  
Carton Weight: 14KG